



OMNİKON

NATIONAL HACKATHON 2026

TEAM NAME

AlphaX duo

MEMBERS

Mounika A , Shanmugapriya S B

THEME

BioTech & HealthTech

PROBLEM STATEMENT

Conventional prosthetics are costly and often lack adaptive grip control, limiting daily functionality for underprivileged amputees. Although affordable and open-source designs exist, they are scattered and hard to discover. Amputees also struggle to find reliable local makers, get proper fitting guidance, and receive aftercare support. As a result, many functional low-cost limbs are either never accessed or abandoned after a short period.

PROBLEM STATEMENT

THE PROBLEM

Choosing a suitable prosthetic is challenging because:

- Users often don't know which prosthetic design suits their needs
- Taking and understanding measurements is difficult without guidance
- Different designs have different mechanical characteristics
- Objective comparison of designs is hard
- Finding trustworthy and verified prosthetic makers is difficult
- Even after selection, simple fitting guidance is limited

Core Problem

- Users need a trusted and understandable pathway from assessment → selection → fitting guidance → reliable makers.

How It Addresses The Problem

- Solves discovery problem → users find suitable low-cost designs easily
- Solves fitting problem → clear guided process reduces pain and abandonment
- Solves trust problem → verified makers + quality signals
- Solves aftercare problem → continuous support and repair guidance

Innovation And Uniqueness

- Mechanical Quality Evaluation of every design (strength, grip usefulness, durability, suitability for Indian activities)
- Trust-first system designed specifically for low-cost market fears
- Fitting and training treated as a core feature, not an afterthought
- Strong emphasis on reducing prosthetic abandonment
- Built with deep understanding of both software and mechanical realities

Introducing ProstheticTrust:

- A digital platform that supports users throughout the prosthetic selection journey.

What we provide:

- Structured Assessment – Collect measurements and user requirements
- Design Recommendation – Match assessment data with design specifications
- Mechanical Quality Score – Clear scores for stability, strength, durability, adjustability & weight
- Fitting Guidance – Simple visual/video-based guidance
- Verified Makers – Discover trusted prosthetic makers
- Feedback & Rating – Build trust through real user experiences

Why ProstheticTrust?

Not just a prosthetic catalogue – it connects the full user journey:

Assess → Understand → Compare → Recommend → Guide → Connect → Feedback

What makes it different?

- Measurement-based decision support
- Mechanical transparency
- Built-in fitting education
- Trust system (Verified makers + Feedback + Ratings)

Features Implemented

Current Working Prototype

Authentication

- Firebase Email/Password Authentication
- Registration, Login, Logout
- User Profile
- Name, Phone, Role, Language
- Stored in Firestore

Assessment

- Affected side, Residual limb length, Forearm & Wrist circumference
- Hand length & width
- Input validation + Scrollable interface

Design Discovery

- Multiple design options
- Mechanical characteristics
- Design selection flow

Navigation Screens

Welcome → Need Selection → Limb Type → Assessment → Review → Recommendation → Guidance



Technical Architecture

Architecture Overview

Flutter App (Dart)

Authentication | Assessment | Review | Recommendations | Fitting Guidance | Maker Discovery | Feedback

Application Logic

- Recommendation Engine | Compatibility Evaluation | Mechanical Scoring | Input Validation
- Firebase Backend
- Firebase Authentication | Cloud Firestore | Secure user-specific data

Technologies

- Flutter
- Dart
- Firebase Authentication
- Cloud Firestore

PROGRESS MODE

Development Progress

Completed

- Flutter project setup
- Firebase integration & Authentication
- User profile + Firestore storage
- Welcome, Need Selection, Limb Type flow
- Assessment screen with validation
- Review Assessment
- Design Recommendation flow
- Fitting Guidance structure

In Development

- Validated recommendation rules
- Mechanical Quality Score
- Measurement guidance videos
- Verified maker system
- Feedback & rating system

Login

Email

Password

Login

Don't have an account? Sign Up

Create Account

Full Name

Email

Phone Number

Password

Role
AMPUTEE

Language
English

Create Account

ProstheticTrust

Test User

Role: AMPUTEE

What would you like to do?

Start Assessment

Take measurements for below-elbow prosthetic

Browse Designs

View available prosthetic designs

Find Makers

Connect with verified makers near you

Assessment

Below-Elbow Prosthetic Assessment

Complete the assessment carefully.

Affected Side

Select side

Residual Limb Length

Length

Forearm Circumference

Circumference

Wrist Circumference

Circumference

Hand Length

Length

Hand Width

HOW TO MEASURE THE BODY FOR FITTING A PROSTHETIC ARM (UPPER LIMB)

Accurate measurement ensures comfort, proper fit and better function.

1 SHOULDER WIDTH



Measure the distance across the shoulders from one edge to the other.

2 CHEST GIRTH



Measure the circumference of the chest, just under the armpits.

3 UPPER ARM LENGTH



Measure from the shoulder point (acromion) to the end of the residual limb.

4 UPPER ARM CIRCUMFERENCE



Measure the circumference of the upper arm at the widest part.

5 ELBOW WIDTH



Measure the width across the elbow (widest part).

6 ELBOW TO END LENGTH



Measure from the elbow joint to the end of the residual limb.

7 FOREARM CIRCUMFERENCE



Measure the circumference of the forearm (widest part).

8 WRIST CIRCUMFERENCE (If Present)



Measure the circumference of the wrist.

9 HAND LENGTH (If Present)



Measure from the wrist crease to the tip of the middle finger.


10 HAND WIDTH (If Present)



Measure the width of the hand across the knuckles.

EXTRA INFORMATION TO COLLECT

- Age
- Weight
- Height
- Cause of limb loss
- Time since amputation
- Activity level & daily needs



GENERAL GUIDELINES

- Use a soft measuring tape.
- Measurements should be taken in centimetres (cm).
- The person should be relaxed and in a natural position.
- Measure the limb (affected side).

STEP 1 — SELECT YOUR AMPUTATION LEVEL

☐ Above the elbow	☐ Below the elbow
☐ At the elbow	☐ At the wrist

STEP 2 — BODY & LIMB MEASUREMENTS

1. Shoulder Width Measure the distance across the shoulders from one prominent shoulder point to the other.	Enter: _____ cm
2. Chest Girth Measure the circumference of the chest, just under the armpits.	Enter: _____ cm
3. Upper Arm Length Measure from the shoulder point (acromion) to the end of the residual limb.	Enter: _____ cm
4. Upper Arm Circumference Measure the circumference of the upper arm at the widest part.	Enter: _____ cm
5. Elbow Width Measure the width across the elbow at the widest part.	Enter: _____ cm
6. Elbow-to-End Length Measure from the elbow joint to the end of the residual limb.	Enter: _____ cm
7. Forearm Circumference Measure the circumference of the forearm at the widest part.	Enter: _____ cm
8. Wrist Circumference If the wrist is present, measure the circumference of the wrist.	Enter: _____ cm
9. Hand Length If the hand is present, measure from the wrist crease to the tip of the middle finger.	Enter: _____ cm

STEP 3 — PHOTOGRAPH THE RESIDUAL LIMB

Upload two clear photographs so the prosthetic specialist can understand the shape of the limb.

FRONT VIEW	SIDE VIEW
Show the entire residual limb from the front. Photo: _____	Show the entire residual limb from the side. Photo: _____

STEP 4 — ADDITIONAL INFORMATION

Age	_____ years	Height	_____ cm
Weight	_____ kg	Affected side	☐ Left ☐ Right ☐ Both
Time since amputation	_____	Activity level	☐ Low ☐ Moderate ☐ High
Cause of amputation	_____		

CURRENT NEEDS — Select all that apply

☐ Daily activities	☐ Work
☐ School / education	☐ Sports / physical activities
☐ Cosmetic appearance	☐ Other: _____

PAIN OR SENSITIVE AREAS

STEP 5 — REVIEW BEFORE SUBMITTING

☒	All measurements are entered in cm.
☒	The tape was not pulled tightly.
☒	Measurements were taken in a relaxed position.
☒	Circumferences and widths were measured at the correct locations.
☒	Front and side photographs are clear and show the residual limb.

SAVE & CONTINUE

CHALLENGES FACED

Challenges & Our Approach

- Measurement Accuracy

Challenge: Users may enter incorrect measurements

Approach: Validation + Instructions + Review screen + Professional verification note

- Recommendation Reliability

Challenge: Arbitrary rules can give poor recommendations

Approach: Validated design specifications + Clear compatibility rules

- Trust

Challenge: Users need confidence in recommendations and makers

Approach: Verified makers + Transparent scoring + Feedback & ratings

- Accessibility

Challenge: Users have different technical abilities

Approach: Simple UI + Visual guidance + Multilingual support + Planned voice assistance

FUTURE ROADMAP

Roadmap

Phase 1 – Functional Prototype (Completed)

Core user journey • Authentication • Assessment • Navigation

Phase 2 – Smart Recommendation (In Progress)

Validated design specifications • Compatibility engine • Explainable mechanical score

Phase 3 – Guidance

Measurement videos • Fitting guidance • Multilingual support • Voice assistance

Phase 4 – Trust Network

Verified makers • Maker profiles • Ratings • User feedback

Phase 5 – Scale

Offline support • Regional maker discovery • Improved accessibility • Data-driven improvements

EXPECTED IMPACT

For Users

- Easier understanding of prosthetic options
- Structured assessment
- Better design comparison
- Accessible educational guidance
- Easier discovery of trusted makers

For Makers

- Structured user requirements
- Better communication
- Opportunity to build trust through verification and feedback

Long-Term Vision

- Build a trusted digital ecosystem connecting users, prosthetic designs, guidance, makers and feedback.